QUICKVUE ADVANCE[®]

G. vaginalis TEST

CLIA Complexity: MODERATE

INTENDED USE

The QuickVue Advance G. vaginalis test is an enzyme activity test for use in the detection of *Gardnerella vaginalis* Proline IminoPeptidase (PIP) activity in vaginal fluid specimens from patients suspected of having bacterial vaginosis. The test is intended for use by healthcare professionals.

SUMMARY AND EXPLANATION

Infectious vaginitis accounts for more than ten million physician office visits each year,¹ with approximately half of all adult women suffering at least one episode. The three most common forms of infectious vaginitis in decreasing incidence are bacterial vaginosis, vulvovaginal candidiasis and trichomoniasis.^{2,3} Because vulvovaginal infections can result in serious clinical sequelae,⁴⁻⁹ vulvovaginal symptoms and signs warrant careful evaluation and appropriate therapy.

The two most common symptoms associated with bacterial vaginosis are vaginal discharge and genital malodor.^{5,10} However, bacterial vaginosis cannot be adequately diagnosed on the basis of symptoms since nearly half of the women who meet the clinical criteria for bacterial vaginosis report no noticeable symptoms.^{4,11}

The Centers for Disease Control and Prevention 1993 Clinical Practice Guidelines state that bacterial vaginosis may be diagnosed by the use of clinical (Amsel) or Gram stain criteria.¹² Culture of *G. vaginalis* is not recommended as a diagnostic method.^{1,8,13,14}

Recent studies have demonstrated that serious clinical sequelae are associated with bacterial vaginosis, including salpingitis,^{4,7} endometritis^{4,7} and post-gynecological surgery infection.^{5,10} In pregnant women, bacterial vaginosis has been associated with placental infection,^{15,16} premature rupture of membranes,^{6,17} pre-term, low birth weight infants,⁹ chorioamnionitis¹⁵ and amniotic fluid infections.¹⁶

The microflora of normal vaginal fluid is characterized by a predominance of aerobic *Lactobacilli*.¹⁸ Bacterial vaginosis results from the replacement of normal vaginal flora with a mixed flora consisting primarily of facultative and anaerobic organisms, including *G. vaginalis*, anaerobic gram negative rods (*Bacteroides spp*, *Prevotella* and *Porphyromonas*), *Mobiluncus spp*. and others.^{6-8,13,19} Of the organisms strongly associated with bacterial vaginosis, only *G. vaginalis* and *Mobiluncus spp*. exhibit PIP activity.^{20,21}

G. vaginalis PIP activity has been colorimetrically detected directly in vaginal fluid specimens from women with bacterial vaginosis using a four hour enzymatic laboratory procedure.²³⁻²⁶ The QuickVue Advance *G. vaginalis* test colorimetrically detects *G. vaginalis* PIP activity directly in vaginal fluid specimens in five minutes or less.

PRINCIPLE OF THE TEST

The QuickVue Advance *G. vaginalis* test system detects PIP activity by using a dried chromogen, Fast Red, to detect hydrolysis of a dried PIP substrate, L-Prolyl- β -naphthylamide. The test employs thin dried films of PIP substrate and chromogen in the test area. The rectangular area immediately adjacent to the test area contains the negative and positive reagent controls. Addition of an undiluted vaginal fluid specimen exhibiting threshold PIP activity hydrolyzes the PIP substrate, resulting in a visible peach-to-pink-to-red color directly on the patient swab within five (5) minutes.

The test detects a threshold activity of thirty (30) picomoles of L-Prolyl- β -naphthylamide hydrolyzed per minute. A positive QuickVue Advance *G. vaginalis* test result indicates that $\geq 2 \times 10^7$ Colony Forming Units (CFU) of *G. vaginalis* are present in the undiluted vaginal fluid.

REAGENTS AND MATERIALS SUPPLIED

- Individually foil wrapped tests (25)
 - ▶ Test contains L-Prolyl- β -naphthylamide and Fast Red.
- Sterile Swabs (25).
- Control Reagent A (7 mL) for use with the controls (1).
- 3-inch swabs for use with the Controls (25).
- Package Insert (1).
- Procedure Card (1).

MATERIALS REQUIRED BUT NOT PROVIDED

- Timer

WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use.
- Do not use kit contents after the expiration date printed on the outside of the kit.
- Dispose of patient samples, containers and used contents in accordance with Federal, State and Local requirements.³²
- Do not remove the test from the foil pouch until ready to use. Once the foil pouch has been opened, the test must be used within 30 minutes. Do not use test if the foil pouch is not sealed.
- Use only the control swabs and control reagent provided in the kit to perform the controls.
- Patient vaginal swabs are not appropriate for any other purpose, including bacterial culture or microscopy, after performing the test.
- Tests with a pink or red color in the test area or control area prior to the application of a sample should be discarded.
- Do not use vaginal fluid specimens that contain blood.
- To obtain accurate results, you must follow the Package Insert instructions.

KIT STORAGE AND STABILITY

Store kit at 2–30°C (36–86°F), out of direct sunlight. Kit contents are stable until the expiration date printed on the outer box carton.

SPECIMEN COLLECTION

Collect an undiluted vaginal fluid sample using **ONLY** a sterile swab provided with the test. The swab can be moistened with vaginal fluid by gently stroking the walls of the vagina. Moisten the swab thoroughly with the sample. The sample should be tested immediately.

QUALITY CONTROL

Built-in Quality Control Features

The QuickVue Advance G. vaginalis test contains built-in positive and negative control features with each test run. For daily quality control, Quidel recommends documenting these controls for the first sample tested each day.

The QuickVue Advance G. vaginalis test contains dried negative and positive controls which monitor test system performance. To perform the controls, add 2–3 drops of the Control Reagent provided in the test kit to the tip of the Control Swab. Firmly rub the Control Swab over the Negative Control Square first. Immediately examine the tip of the swab for the presence or absence of a pink-to-red color. Next, firmly rub the same Control Swab over the Positive Control Square. Immediately examine the tip of the swab for the presence or absence of a pink-to-red color.

Negative Reagent Control

The dried Negative Reagent Control area contains Fast Red. It should not produce any pink-to-red color on the Control Swab. If a pink-to-red color forms on the Control Swab, this indicates a potential failure of the test system. Test results should be disregarded and/or repeated.

Positive Reagent Control

The dried Positive Reagent Control area contains a reagent (3-amino-2-naphthoic acid) that reacts chemically with the Fast Red chromogen to produce a visible pink-to-red color directly on the Control Swab. It should always produce a pink-to-red color on the Control Swab. If no pink-to-red color forms on the Control Swab, this indicates a potential failure of the test system. Test results should be disregarded and/or repeated.

NOTE: The formation of a peach-to-pink-to-red color on the patient swab is qualitative and should not be compared in terms of color hue to intensity to color formed on the positive reagent control swab.

TEST PROCEDURE

Note: The test result is on the swab. Do not discard swab until results are read.

1. Remove test from the foil pouch just prior to use.
2. Using a circular motion firmly rub the moistened patient swab over the entire surface of the test area approximately ten times.
3. Examine the tip of the swab for a peach-to-pink-to-red color. Lay the patient swab on the QuickVue Advance G. vaginalis test with the tip over the test area.
4. Read result within five (5) minutes.

INTERPRETATION OF RESULTS

Results are determined by the presence or absence of color on the patient swab.

Positive Result

The **presence of a visible peach-to-pink-to-red color on the tip of the patient swab** within five (5) minutes is a positive result. A positive test result indicates the detection of PIP activity in the vaginal fluid specimen associated with *Gardnerella vaginalis* and is suggestive of bacterial vaginosis if consistent with other clinical signs and symptoms, including elevated pH; an amine odor; thin, homogeneous, adherent vaginal discharge; and clue cells.

Negative Result

The **absence of a visible peach-to-pink-to-red color on the tip of the patient swab** at five (5) minutes is a negative result. A negative test result indicates that PIP activity associated with the presence of *Gardnerella vaginalis* is absent or below the detectable level in the sample collected. Gram stain criteria as well as clinical signs and symptoms should be considered before suggesting that the woman does not have bacterial vaginosis.

LIMITATIONS

- Test results may be affected by improper specimen collection and/or handling. A negative test result does not exclude the possibility of bacterial vaginosis.
- Mixed infections may occur. Therefore, a test indicating the presence of *G. vaginalis* does not rule out the presence of yeast, *Trichomonas vaginalis* or other organisms.
- The presence or absence of *G. vaginalis*, although suggestive, is not diagnostic for bacterial vaginosis. Results should be interpreted in conjunction with other clinical data including vaginal discharge characteristics, pH, amine odor, and clue cells.

EXPECTED VALUES

G. vaginalis can be detected in approximately fifty percent of normal women.¹³ However, the vaginal fluid Geometric Mean Concentration (GMC) of *G. vaginalis* is more than 100 times higher in women with bacterial vaginosis than in normal women (1.3×10^3 versus 9.7×10^5 CFU/ml).¹³ This concentration of vaginal *G. vaginalis* has been associated with bacterial vaginosis in a study employing quantitative *G. vaginalis* culture and a research oligonucleotide probe.²²

PERFORMANCE CHARACTERISTICS

The performance of the QuickVue Advance *G. vaginalis* test to detect *G. vaginalis* PIP activity in vaginal fluid specimens collected from symptomatic and asymptomatic patients, was determined in a multi-center study conducted at three geographically separate clinical sites. A total of 640 women were enrolled and consecutively tested. Study site personnel performed the following: the QuickVue Advance *G. vaginalis* test; the four clinical criteria - vaginal fluid homogeneity, pH test, whiff test and clue cells on microscopy;²⁷ wet mount microscopy for candidiasis and trichomoniasis; and prepared a microscope slide for subsequent Gram stain analysis for testing at a CLIA certified laboratory. The Gram stain score was classified as follows: normal (0–3), intermediate (4–6), and bacterial vaginosis positive (7–10). The *G. vaginalis* morphotype Gram score (0–4+) was also determined. *G. vaginalis* morphotypes are scored from 0 (no morphotypes) to 4+ (greater than 30 morphotypes).¹⁴

The results of the QuickVue Advance *G. vaginalis* test were compared to the results of the Gram stain for diagnosing bacterial vaginosis. Discordant results were reconciled with the results of the clinical or Amsel criteria. Of the 451 symptomatic women, only 62% were diagnosed as having at least one form of infectious vaginitis (bacterial vaginosis, candidiasis and/or trichomoniasis). Of the 185 asymptomatic women, 25% were diagnosed as having at least one form of infectious vaginitis. The sensitivity (95.1%) and specificity (98.6%) of the QuickVue Advance *G. vaginalis* test in the asymptomatic population did not differ to a statistically significant extent from the sensitivity (90.9%) and specificity (97.1%) in the symptomatic population.

Comparison to Gram stain and the Clinical or Amsel Criteria

The performance of the QuickVue Advance *G. vaginalis* test was compared to the two reference methods for diagnosing bacterial vaginosis (Gram stain and the clinical or Amsel criteria). Only those women with at least three of the four clinical criteria or a Gram stain score of 7 or above were considered positive for bacterial vaginosis (BV). All Gram stain intermediate (Gram score of 4–6) test results were treated as bacterial vaginosis negative.

Out of the 637 specimens tested, the two reference methods (Gram stain and clinical criteria) agreed in 552 cases (86.7%). In the 222 specimens where both reference methods produced positive tests for bacterial vaginosis, the QuickVue Advance *G. vaginalis* test was positive in 201 cases (90.5%). Similarly, in the 330 specimens where both reference methods produced negative test results for bacterial vaginosis, the QuickVue Advance *G. vaginalis* test was negative in 321 cases (97.3%). Overall agreement between the QuickVue Advance test and the results of the combined Gram stain and clinical criteria was 96.2%.

TABLE 1

Comparison of the QuickVue Advance G. vaginalis test vs. Gram Stain and Clinical Criteria

Reference Method		QuickVue Advance Test Result	Number of Patients
Gram Stain	Clinical (Amsel) Criteria		
+	+	+	201
+	+	—	21
+	—	+	22
+	—	—	36
—	+	+	7
—	+	—	20
—	—	+	9
—	—	—	321
Total			637

Comparison To Gram Stain Results

QuickVue Advance G. vaginalis test results were compared to the Gram stain results for diagnosing bacterial vaginosis (**Table 2**). Discordant results were reconciled using the clinical or Amsel criteria. Overall agreement was 93.9%.

TABLE 2

QuickVue Advance G. vaginalis test vs. Amsel-Reconciled Gram Stain

Comparative Method	Result	Reference Method Amsel-Reconciled Gram		N
		Positive	Negative	
QuickVue Advance Test Result	Positive	230	9	637
	Negative	21	377	

	Analytical %	Standard Deviation %	Lower Bound 95% One Sided CI
Sensitivity	91.6	1.7	88.8
Specificity	97.7	0.8	96.4
Pos Predictive Value	96.2	1.2	94.2
Neg Predictive Value	94.7	1.1	92.9
Overall Agreement	95.3	0.8	93.9

Comparison To Clinical Criteria Results

QuickVue Advance G. vaginalis test results were also compared to the results of the clinical criteria to diagnose bacterial vaginosis. Discordant results were reconciled using Gram stain. When the performance of the QuickVue Advance G. vaginalis test was compared with the Gram-reconciled clinical criteria, identical performance characteristics were found to that seen with the Amsel-reconciled Gram stain. The QuickVue Advance G. vaginalis test produced the following results: sensitivity = 91.6%; specificity = 97.7%; positive predictive value = 96.2%; negative predictive value = 94.7%; and overall agreement = 95.3%. The standard deviations and the lower bound of the 95% one sided confidence intervals were also identical to those seen in Table 2.

Comparison To Gram Stain Morphotypes

The positivity of the QuickVue Advance G. vaginalis test with specific vaginal fluid levels of Gram stain morphotypes determined by microscopy was established (**Table 3**).

TABLE 3

Positivity of Three Test Methods with the Gram Stain Morphotype Score

Reference Method	Number of Specimens	Gram Stain Morphotype Score				
		0+ n=248	1+ n=58	2+ n=37	3+ n=14	4+ n=279
Gram Stain BV	280/636 44.0%	0/280 0%	0/280 0%	0/280 0%	1/280 0.4%	279/280 99.6%
Clinical Criteria	249/636 39.2%	7/249 2.8%	10/249 4.0%	7/249 2.8%	4/249 1.6%	221/249 88.8%
QuickVue Positive	239/636 37.6%	3/239 1.3%	3/239 1.3%	8/239 3.3%	3/239 1.3%	222/239 92.9%

Gram Stain Of 636 vaginal fluid specimens analyzed, 280 (44%) had a Gram stain score of 7–10. Of these, 279 had a 4+ *G. vaginalis* morphotype score (99.6%). The remaining Gram stain positive specimen had a *G. vaginalis* morphotype score of 3+. Of the 357 specimens having a *G. vaginalis* score of 0–3+, only one (0.3%) had a Gram stain score of 7–10.

Clinical Criteria Of the 636 vaginal fluid specimens analyzed, 249 (39.2%) had three of the four criteria. Of these, 221 had a *G. vaginalis* morphotype score of 4+ (88.8%). Of the 357 specimens having a *G. vaginalis* morphotype score of 0-3+, only 28 (7.8%) had three of the four criteria.

QuickVue Advance G. vaginalis test Of the 636 vaginal fluid specimens analyzed, 239 (37.6%) were QuickVue Advance G. vaginalis test positive. Of these, 222 had a *G. vaginalis* morphotype score of 4+ (92.9%). Of the 357 specimens having a *G. vaginalis* morphotype score of 0-3+, only 17 (4.8%) were QuickVue Advance test positive.

Mixed Clinical Diagnosis

Bacterial vaginosis is often asymptomatic; and even when symptoms are present, they correlate poorly with a clinical diagnosis of vaginitis/vaginosis.^{4,11,25} The relative lack of specificity of symptoms precludes a differential diagnosis based on symptoms. In the clinical study, bacterial vaginosis (BV), candidiasis (CD) and trichomoniasis (TR) were diagnosed on the basis of vaginal discharge, pH, amine odor and wet mount results.

Table 4 shows the results of different test methods in patients stratified by clinical diagnosis.

TABLE 4

Test Results In Patients Stratified By Clinical Diagnoses

No. of Patients	Test Result				
	Amsel Criteria BV	Wet Mount		Gram Score 7-10	QuickVue Advance
		CD	TR		
0	+	+	+	0	0
9	+	+	—	8	7
23	+	—	+	18	16
2	—	+	+	1	1
11	—	—	+	6	3
66	—	+	—	9	4
217	+	—	—	196	185
308	—	—	—	41	23
636				279	239

Quantitative PIP Specific Enzyme Activity in Clinically Characterized Vaginal Fluid Specimens: To verify the utility of vaginal fluid PIP activity as an aid in the diagnosis of BV, the PIP activity in 250 vaginal fluid samples was quantitated. Results of these analyses, and the correlation between median vaginal fluid PIP activity and four reference criteria (Gram stain, Amsel criteria, QuickVue Advance G. vaginalis test and Amsel-reconciled Gram Stain) are found in **Table 5**.

TABLE 5

Comparison of Quantitative PIP Activity in Vaginal Fluid Samples with Reference Methods

Reference Method	Results	No. Tested	Median PIP SEA	Ratio SEA Positive	
				SEA Positive	SEA Negative
Gram Stain	Positive	108	5.44	34	
	Negative	142	0.16		
Amsel Criteria	Positive	106	5.18	30	
	Negative	144	0.17		
QuickVue Advance Test Results	Positive	99	5.49	34	
	Negative	151	0.16		
Amsel-Reconciled Gram Stain	Positive	103	5.49	34	
	Negative	147	0.16		

SEA = Specific Enzyme Activity
(nmoles L-Prolyl-β-naphthamide hydrolyzed per minute per mL of vaginal fluid)

The median PIP specific enzyme activity (SEA) of 108 Gram stain positive specimens was 34 fold higher than the 142 negative specimens. Similarly, the median PIP SEA of 106 Amsel criteria positive specimens was 30 fold higher than the 144 negative specimens. The median PIP SEA of the 103 Amsel-reconciled Gram stain positive specimens was also 34 times higher than that of the 147 negative specimens. Specimens were also grouped on the basis of their QuickVue Advance test results (99 positive and 151 negative). The ratio of PIP median SEA in QuickVue Advance positive to negative specimens was 34. Using Mann-Whitney Two-Sample Rank Test statistics, the differences in median PIP values between the positive and negative populations, regardless of the criteria used, were found to be significantly different (P = < 0.0001).

NON-CLINICAL RESULTS

Specificity

A total of 32 culture grown microbial isolates typically found in vaginal fluid from normal women and women with vaginal and cervical infections were used to challenge the microbial specificity of the QuickVue Advance test (**Table 6**). All cultured organisms were tested at undiluted cell densities of at least 3 x 10⁷ Colony Forming Units (CFU/mL). Where practicable, microbial concentrations of up to 6 x 10⁹ CFU/mL were tested. The challenge cell concentrations, therefore, were 7.5–1500 times higher than that required to produce a positive QuickVue Advance test result with *G. vaginalis* cells. All of the organisms tested except *Gardnerella vaginalis*, *Mobiluncus curtisii* and *Mobiluncus mulieris* were negative in the QuickVue Advance test. Both *G. vaginalis* and *Mobiluncus spp.* are associated with bacterial vaginosis and exhibit PIP activity.^{20,21,23–26} The concentration of undiluted

Mobiluncus spp. required to produce a positive QuickVue Advance test result was $> 4 \times 10^8$ CFU/mL or 20 fold higher than that required to produce a positive QuickVue Advance test result with *G. vaginalis* cells.

TABLE 6

Organisms Tested in Specificity Challenge

<i>Acinetobacter</i> spp.	<i>Aclinomyces</i> spp
<i>Bacteroides levii</i>	<i>Bacteroides ureolyticus</i>
<i>Campylobacter</i> spp.	<i>Candida albicans</i>
<i>Chlamydia trachomatis</i>	<i>Clostridium</i> spp.
<i>Corynebacterium</i> spp.	<i>Cryptococcus</i> spp.
Enterobacteriaceae	<i>Enterococcus</i> spp.
<i>Fusobacterium nucleatum</i>	<i>Gardnerella vaginalis</i>
<i>Haemophilus ducreyi</i>	<i>Lactobacillus acidophilus</i>
<i>Mobiluncus curtissii</i>	<i>Listeria</i> spp.
<i>Mobiluncus hominus</i>	<i>Mycobacterium</i> spp.
<i>Peptostreptococcus telradius</i>	<i>Neisseria gonorrhoeae</i>
<i>Prevotella bivia</i>	<i>Peptostreptococcus anaerobius</i>
<i>Propionibacterium</i> spp.	<i>Prevotella disians</i>
<i>Staphylococcus</i> spp.	<i>Pseudomonas</i> spp.
<i>Trichomonas vaginalis</i>	<i>Streptococcus</i> spp.

Analytical Sensitivity

The QuickVue Advance *G. vaginalis* test can detect $\geq 2 \times 10^7$ CFU/mL of *G. vaginalis* in the undiluted vaginal fluid. The test detects a threshold activity of thirty (30) picomoles of L-prolyl- β -naphthylamide hydrolyzed per minute.

Reproducibility

The reproducibility of the QuickVue Advance *G. vaginalis* test was evaluated at three sites by three different types of users: 1) at 3 physician offices by 4 first time users; 2) at the 3 clinical sites by 11 of the study participants; and 3) in-house by 3 trained users. The samples tested contained different levels of PIP activity corresponding to a negative test result and a positive test result. A total of 119 observations were made by the first time users. They interpreted 30/30 (100%) of all negative samples correctly and 87/89 (98%) of all positive samples correctly. A total of 94 observations were made by the clinical site personnel. They interpreted 30/30 (100%) of all negative samples correctly and 62/64 (97%) of all positive samples correctly. A total of 240 observations were made by the trained users. They interpreted 60/60 (100%) of all negative samples correctly and 180/180 (100%) of all positive samples correctly.

The results of these reproducibility studies confirm that there is virtually no difference between first time users, trained clinicians and trained in-house users to perform and interpret the QuickVue Advance *G. vaginalis* test.

ASSISTANCE

If you have any questions regarding the use of this product, please call Quidel's Technical Support number, (800) 874-1517 (toll-free in the U.S.A.) or (858) 552-1100, Monday through Friday, between 7:00 a.m. and 5:00 p.m. Pacific Time, U.S.A. If outside the United States, contact your local distributor or technicalsupport@quidel.com.

REFERENCES

1. Sobel, J. D. 1990. *Obstet. Gynecol. Clin. N.A.*, 17: 851-879.
2. Holmes, K.K. In: Holmes, K.K., Mardh, P-A., Sparling, P.F., et al eds. *Sexually Transmitted Diseases*, McGraw-Hill, New York. 1990. Chapter 46; 527-545.
3. Wasserheit, J.N. and Holmes, K.K. 1985. *Em. Med. Clin. North America*, 3: 47-74.
4. Eschenbach, D.A., Hillier, S., Critchlow, C., Stevens, C., DeRouen, T., Holmes, K.K. 1988. *Am. J. Obstet. Gynecol.*, 158 (4): 819-828.
5. Soper, D.E., Bump, R.C., Hurt, W.G. 1990. *Am. J. Obstet. Gynecol.*, 163: 1016-1023.
6. Martius, J., Krohn, M., Hillier, S., Stamm, W., Holmes, K.K. and Eschenbach, D.A. 1988. *Obstet. Gynecol.*, 71: 89-95.
7. Holst, E. 1990. *J. Clin. Microbiol.*, 28: 2035-2039.
8. Spiegel, C.A., Ansel, R., Holmes, K.K. 1983. *J. Clin. Microbiol.*, 18: 170-177.
9. Hillier, S.L., Nugent, R.P., Eschenbach, D.A., Krohn, M.A., Gibbs, R.S., Martin, D.H., Cotch, M.F., Edelman, R., Pastorek, J.G., Rao, A.V., McNellis, D., Regan, J.A., Carey, C., Klebanoff, M.A. 1995. *N. Eng. J. Med.*, 333: 1737-1742.
10. Larsson, P.G., Platz-Christensen, I.J., Forsum, U., Pahlson, C. 1991. *Obstet. Gynecol.*, 77.
11. Biswas, M.K. 1993. *Clin. Obstet. Gynecol.*, 36: 166-176.
12. Centers for Disease Control. *Sexually Transmitted Diseases Clinical Practice Guidelines*. US Department of Health and Human Services. 1991: 111-3.
13. Hillier, S.L. and Holmes, K.K. In: Holmes, K.K., Mardh, P-A., Sparling, P.F. et al. eds. *Sexually Transmitted Diseases*. McGraw-Hill, New York. 1990. Chapter 47; 547-559.
14. Nugent, R.P., Krohn, M.A., Hillier, S.L. 1991. *J. Clin. Microbiol.*, 29: 297-301.
15. Hillier, S., Martius, J., Krohn, M., Kiviat, N., Holmes, K.K. and Eschenbach, D.A. 1988. *N. Engl. J. Med.*, 319: 972-978.
16. Gravett, M.G., Nelson, H.P., DeRouen, T., Critchlow, C., Eschenbach, D.A., Holmes, K.K. 1986. *JAMA*, 256: 1899-1903.
17. Martius, J., Eschenbach, D.A. 1990. *Arch Gynecol. Obstet.*, 247: 1-13.
18. Spiegel, C.A., Ansel, R., Eschenbach, D., Schoenknecht, F. and Holmes, K.K. 1980. *N. Engl. J. Med.*, 303: 601-607.
19. Krohn, M.A., Hillier, S.L., Eschenbach, D.A. 1989. *J. Clin. Microbiol.*, 27 (6): 1266-71.

-
20. Balows, A., Hausler, Jr., W.J., Herrmann, K.L., Isenberg, H.D., Shadomy, H.J., (Eds.) Manual of Clinical Microbiology, American Society for Microbiology, Washington, DC 1991 Fifth Edition.
 21. Schoonmaker, J.N., Lunt, B.D., Lawelin, D.W., French, J.I., Hillier, S.L. and McGregor, J.A. 1991. Am.J. Obstet. Gynecol., 165: 737-742.
 22. Sheiness, D., Dix, K., Watanabe, S., Hillier, S.L. 1992. J. Clin. Microbiol., 30: 642- 648.
 23. Thomason, J.L., Gelbart, S.M., Wilcoski, L.M., Peterson, A.K., Jilly, B.J. and Hamilton, P.R. 1968. Obstet. Gynecol., 71: 607-811.
 24. Thomason, J.L., Gelbart, S.M., Anderson, R.J., Walt, A. K., Osypowski, P. J. and Broekhuizen, F.F. 1990. Am. J. Obstet. Gynecol., 162: 155-160.
 25. Livingood, C. H., Thomason, J.L. and Hill, G.B. 1990. Am. J. Obstet. Gynecol., 163: 515-20.
 26. James, J.A., Thomason, J.R., Gilbert, S.M. 1992. Am. J. Obstet. Gynecol., 166; 859-863.
 27. Amsel, R., Totten, P.A., Spiegel, C.A., Chen, K.C.S., Eschenbach, D.A., Holmes, K.K. 1983. Am.J. Med. 74: 14-22.
 28. Rotimi, V.O., Yakubu, Z., Abudu, O.O., Banjo, T.O. 1991. J. Med. Microbiol., 35: 103-6.
 29. Hillier, S.L., Krohn, M.A., Nugent, R.P., Gibbs, R.S. 1992. Am. J. Obstet. Gynecol., 166: 938-44.
 30. Mazzulli. T., Simor, A.E., Low, D.E. 1990. J. Clin. Microbiol., 28: 1,506-8.
 31. Joesoef, M.R., Hillier, S.L., Suharo, J. and Linnam, M. 1991. J. Clin. Microbiol., 29:1730-1.
 32. Recommendations for the Prevention of HIV Transmission in Health Care Settings, Morbidity and Mortality Weekly Report, Centers for Disease Control and Prevention, August 21, 1987.

Covered by U.S. Patent No. 5,571,684; European Patent No. 0 791 071; other patents pending.

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Quidel Corporation
Worldwide Headquarters
10165 McKellar Court
San Diego, CA 92121 USA
www.quidel.com

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