

QUICKVUE®

Influenza test



Test and Treat Today!

Diagnosing influenza can be confusing. One third of all upper respiratory infections are caused by influenza. This leads to over 36,000 deaths, 114,000 hospitalizations and accounts for over \$12 billion spent in lost work and related medical expenses each year in the U.S.

Quidel's CLIA-waived QuickVue® Influenza test provides you with test

results that allow you to Test and Treat in less than 10 minutes.

Prompt diagnosis and treatment of influenza is a solution to the overuse of antibiotics, allows for timely use of antivirals, and prevents healthcare costs related to unnecessary testing and treatment.



● RAPID DIAGNOSIS AND TREATMENT

- Results in 10 minutes or less (right in the exam room) for better patient management decisions.
- Most patients tested can be diagnosed and treated in just one office visit.

● EASY TO USE

- 1 step (nasal wash/aspirate), 3 steps (nasal swab), 1 reagent, 30 seconds hands-on time.
- First CLIA-waived flu test.

● ACCURATE DETECTION

- QuickVue® Influenza test detects both types A and B viral antigens.
- High accuracy.

● CONVENIENT PACKAGING



- Each kit contains everything needed to perform the test, including both positive and negative controls.
- Available in two kit sizes: 10-test (individual trays) and 25-test (bulk).
- Room temperature storage with up to 24-month shelf life.
- Transport swab set also available separately.

For more information on QuickVue® Influenza test kits, visit our website at www.quidel.com, or call your Quidel Account Manager at 1.800.874.1517

QUIDEL®
CORPORATION

Helping Women and Their Families Live Healthy Lives

Nasal swab procedure:





Sensitivity: 73% (95% C.I. 67%-78%)
Specificity: 96% (95% C.I. 93%-98%)

Pred. Value (+): 92%
Pred. Value (-): 85%

Accuracy: 87%

Nasal wash/nasal aspirate procedure:



Sensitivity: 81% (95% C.I. 72%-88%)
Specificity: 99% (95% C.I. 93%-99%)

Pred. Value (+): 96%
Pred. Value (-): 93%

Accuracy: 94%

To order: Call your Quidel Account Manager at 1.800.874.1517

10 Test Kit: Catalog #00324
25 Test Kit: Catalog #00317
Transport Swab Set: Catalog #20171

QUIDEL[®]
CORPORATION

Helping Women and Their Families Live Healthy Lives

Quidel Corporation • Worldwide Headquarters • 10165 McKellar Court • San Diego, CA 92121 USA • www.quidel.com
For more information on our products and other Quidel locations, visit our web site.

Distributor Information: