

Why More Professionals Are Choosing INRatio...

The INRatio PT/INR testing system combines practical convenience with reliable results making it the perfect in-office solution for anticoagulation patient management.

On Board Quality Control...

Rigorous quality control checks are performed automatically with every test strip, ensuring a reliable result.



No Refrigeration Required...

Save space! No waiting for strips to warm up. Store at room temperature for up to 12 months!



Small Sample Size...

1 drop of blood (15µL) applied directly to the strip performs both the test and quality controls.



Reimbursement Information*

CPT Code	Description
85610QW	Prothrombin Time (The QW modifier indicates a CLIA waived test)
99201-99205	E&M of NEW Patients
99211-99215	E&M of ESTABLISHED Patients

NOTE: Reimbursement may vary depending on your location/plan. It is the responsibility of the health care provider to properly code and seek reimbursement.

References

¹ The Seventh ACCP Conference on Anticoagulation Therapy: Evidence-Based Guidelines, Supplement to Chest, Volume 126, No. 3 (Supplement), September 2004: 2085-2095.

² Chiquette E., Amato M., Bussey H., Comparison of an Anticoagulation Clinic with Usual Medical Care. Arch Intern Med. 1998; 158: 1611-1647.


About HemoSense

HemoSense was founded to develop next generation point-of-care blood coagulation testing devices for the hemostasis markets. We seek to empower healthcare professionals and patients with the confidence, control, and convenience needed to safely and effectively manage oral anticoagulation therapy.

HemoSense®...

HemoSense, Inc.
651 River Oaks Parkway
San Jose, CA 95134
Phone 408-719-1393
Fax 408-719-1184
moreinfo@hemosense.com
www.hemosense.com

PN 5500279 RevB ©HemoSense. All rights reserved.



CLIA Waived

Actual Size



Real time, On Board Quality Controls. Every time you test.

Quality control is essential to PT/INR testing. With INRatio, 2 levels of QC are performed automatically every time you test. **No more mixing separate liquid controls or using additional test strips for QC testing. That saves time and money.** For added confidence, the QC results appear on screen.

Total System Confidence...

INRatio also:

- Performs electronic self checks.
- Ensures test strip is inserted properly.
- Verifies operating temperature.
- Checks for sufficient sample size.

With all these integrated safety checks, you can be confident that your patients are being managed with reliable results every time you test.

Each test strip performs a normal (QC1) and a therapeutic (QC2) check on every patient sample.

Easy Test Procedure...

Simply apply 1 drop of blood (15µL) directly to the strip to perform both the test and quality controls. No cups or cuvettes needed.



Accurate Results...

The INR is calculated with a low ISI, meeting the recommendations of the American College of Chest Physicians, College of American Pathologists, and the World Health Organization.¹

INRatio Supports Your Goals of Improving Patient Care and Practice Efficiencies.

INRatio will enable you to minimize the challenges of managing your oral anticoagulation patients.

Improve patient satisfaction...

Send patients home instead of to the lab for a PT/INR test.

Eliminate unnecessary delays...

Get results in less than 2 minutes instead of hours or even days.

Increase compliance...

Counsel your patients face-to-face with dosage changes for improved TTR (Time in Therapeutic Range).²



Maximize practice efficiencies...

Reduce time-consuming "call-backs", freeing your staff for more productive duties.

E-Trainer...

A self-paced, interactive, web or CD based educational tool. Learn more at www.hemosense.com today!

INRatio Pro Kit Cat #0100004

- 1 INRatio Monitor
- 1 Power Supply
- 50 Blood Lancets
- 1 User's Guide
- 1 Quick Reference Guide
- 1 Training Video

Test Strips Cat #0100139
48 Strip Box